

Detection and confirmation of Androgen receptor -ARV7


Gold Standard Accuracy — Verifiable Results Format — Multiplexed — High Throughput-Rapid TAT — Cost-Effective

Testing for Abiraterone (Zytiga®) resistance prostate cancer

ARV7- Tumorplex™ (For solid tumors and plasma)

Key Features

- Detection and confirmation Arv-7 variant.
- Highest sensitivity **Mutant: Wild >1:100,000 vs 1:1000**
- Wild type as internal control without suppression.
- No need to **enrich tumor cells** = cost savings.
- Small sample volume.
- **LOD: Solid tumors 15 copies of genomic DNA
Plasma: Copies of 110 base RNA**



wiseGEEK

