


Detection and Confirmation of EGFR del 747_A750

Gold Standard Accuracy Verifiable Results Format High Throughput-Rapid TAT Cost-Effective


EGFR del- Tumorplex™ (for solid tumors and plasma)

KEY FEATURES

- Detection and confirmation EGFR del 747_A750.
- Wild type as internal control without suppression.
- No need to enrich tumor cells = cost savings.
- Small sample volume.
- LOD: Solid tumors 1 copy



Indicator – Wild Type



Indicator – Wild Type

Indicator 747_A750 deletion

